Scintillomètre SPP2 NF

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1 - USE :

The SPP 2 NF Portable Prospecting Scintillometer is designed to detect and measure the gamma radiation given off by radioactive substances. It represents a new version of the SPP 2 SRAT and is produced throughout in printed circuits of the removable and plugin type, which greatly facilitates repair and improves its strength. Its electronic assembly, which is schematically identical with that of the SPP 2, is such that it can be included in the precision equipment category. The SPP 2 NF also contains a sound alarm device with a threshold and frequency that can be varied according to the strength of the radiation.

2 - PRINCIPAL CHARACTERISTICS :

The SPP 2 NF Scintillometer basically comprises a scintillator (Na I.Tl) photomultiplier assembly used as a detector, an impedance matcher, an electronic flip flop, a direct-reading integrator, a stabilised transistor power-supply and a pulse-generating circuit the sound alarm.

Power Supply

The SPP 2 NF is powered by three 1,5 V torch batteries arranged in series. (R20)

Operating life :

Its operating life is 80 hours in intermittent operation (20 times 4 hours). Its life is 40 hours in continuous operation with the sound alarm disconnected (threshold button hard to the right) and about 20 hours with the alarm operating, using R20 batteries. These times can be increased to 50 and 25 hours respectively with R 20S batteries.

Weight :

Its weight is 3.600 kgs in running order, with batteries.

Dimensions :

The external dimensions are :

Length	32 cm
Width	13 cm
Height	21 cm

Water tightness :

The apparatus is designed for an immersion under 1 m of water.

Accuracy :

- In normal conditions of measurement, that is to say : Temperature $\pm 20^{\circ}\text{C} \pm 10^{\circ}\text{C}$ and gamma radiation with an energy above 30 keV : $\pm 10\%$.*
- Statistical error products by the random character of the information: see fluctuations to paragraph «Sensitivity» (below).*

Variation products by temperature influence :

In normal range of utilization, that is to say, - 20°C to $+ 60^{\circ}\text{C}$: $\leq \pm 20\%$.

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Accessories :

A «Pulses» output ($\square$) permits a headphone branching. A variable direct current output ($=$) on a 1 000 Ω impedance permits to record the different ranges, by means of a potentiometer recorder ETP 3 type S.R.A.T. for instance.

Sensitivity :

(altered by changing the position of the switch lever)

To each sensitivity there corresponds an integration constant value giving a specific rate of fluctuation at 1/3 of the scale, as indicated below.

Sensitivity in real hits per second	Integration Constant	Fluctuation at 1/3 of scale
1) 0 - 15 000	0.06 s	4.1 %
2) 0 - 5 000	0.06 s	5.9 %
3) 0 - 1 500	0.4 s	5 %
4) 0 - 500	0.6 s	7 %
5) 0 - 150 fast	0.6 s	13 %
6) 0 - 150 slow	4 s	5 %

Another position of the switch marked «Contrôle» enables the conditions of the batteries to be checked before radiation measurements are performed. A set of batteries is unusable when, in the «Contrôle» position, the needle of the measuring device does not reach the blue area.

Another last, position of the switch marked «Arrêt» causes the power supply battery circuit to open and short circuits the measuring device.

Built-in sound alarm :

Sound-alarm threshold adjustable by a milled external button between one-twentieth and maximum of the scale on each of the instrument's sensitivities - Starting from the threshold selected, the frequency of the sound emitted is largely proportional to the intensity of the incident radiation.

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1 - APPEARANCE :

1-1 - EXTERNAL LAYOUT :

The SPP 2 NF scintillometer comes in the form of a duralinox box coupled to a case containing the batteries and the sound alarm. This assembly supports a pistol-shaped block fitted with a grip. The pistol is electrically connected to the case by a cable measuring about 1,5 metres.

At the side of the box, attachments and passages for straps are provided for wearing the equipment at the waist.

On the box there are the sensitivity-switch control button and the outlet points for the pulses and for the integrated continuous voltage.

On the back of the battery container and to the right there are the sound alarm adjustment button and, on the left, the sound outlet aperture.

1-2 - INTERNAL LAYOUT :

a) Box :

The whole electronic part is attached to one of the sides of the box. The electronic assembly is composed of five sub-assemblies :

- 1) *An H.V (1 000 v + 500 v) and M.V. (12 v) power supply sub assembly containing three printed circuit mounting-plates attached to one another by angle pieces and electrically linked by wires. One of the three plates carries the set of connecting plugs.*
- 2) *An electronic-rocker sub-assembly on a single plate, together with its set of connecting plugs. It is linked by three angle-pieces to the power-supply sub-assembly but is electrically independent*
- 3) *A buzzer sub-assembly which can be plugged in separately and which controls the sound alarm.*
- 4) *An integration sub-assembly, which can also be plugged in separately*
- 5) *An inter connection sub-assembly containing a double-sided printed-circuit plate and bearing the sensitivity switch, into which sub-assembly the four previous ones are plugged.*

By means of a torus of flexible wires leading to a set of plugs on the base of the box, the electronics can be easily extracted from the box and used outside it.

b) Pistol body (Probe) :

This comprises :

- 1) *The photomultiplier-scintillator coupling with dynode power supply.*
- 2) *The impedance-matching sub-assembly, in printed circuit form.*
- 3) *The measuring device.*

Note : The position of the scintillator crystal in the probe is marked on the outside by a red circle and ring.

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2 - SCHEME AND OPERATION (see general diagram Appendix 1)

2-1 - DETECTOR CIRCUIT :

The detector circuit of the SPP 2 NF comprises :

- a) A 1 x 1.5 inch (SC 1) NaI TI (sodium iodide activated with thallium) scintillator crystal.*
- b) A XP2008 photomultiplier powered by the H.V. supply described later. (PM1)*

An incident particle of energy E_1 which penetrates into the scintillator will dissipate in the latter an energy E_2 less than or equal to E_1 . This absorbed energy E_2 will be restored by the scintillator in the form of N_0 photons of energy E such that $N_0.E$ is less than or equal to E_2 .

Of these N_0 photons produced in the scintillator, a number N_1 less than or equal to N_0 arrives on the photocathode of the photomultiplier and proceeds to extract N_2 electrons from it by photoelectric effect. Of these N_2 electrons, N_3 less than or equal to N_2 are collected by the first dynode of the photomultiplier and are then multiplied up to the anode, where N electrons are found. This charge flows into the anode circuit and gives a voltage pulse collected on resistance R_{13} .

This voltage pulse has a level directly proportional to the energy lost in the crystal by the incident particle. In the case of counting pulses of very different levels are therefore obtained at the terminals of resistance R_{13} whence the necessity to calibrate them in amplitude and time before integrating them, this is the function of the counting circuit.

In order to transmit correctly the pulses from the photomultiplier to the counting circuit via a cable, it is necessary to lower the impedance at the cable entrance.

2-2 - IMPEDANCE MATCHING CIRCUIT :

This circuit is established by means of an NPN (Q1) transistor mounted as a common collector.

The positive pulses collected at the terminals of R_{13} are applied to its base and collected on the slider of a potentiometer R_{30} connected in its emitter (R_{30} is situated on the buzzer plate).

They are then transmitted across the linking condenser C_{15} to the counting circuit (flip flop plate).

2-3 - COUNTING CIRCUIT :

The basic component of the circuit for counting the pulses from the detector circuit is a quick return, monostable electronic flip flop.

The function of this rocker is to send constant amplitude signals into the integrating circuit, starting from short pulses originating in the photomultiplier. These calibrated signals are integrated in the reading circuit described in the next paragraph.

The rocker contains two PNP transistors, Q 3 and Q 4. These two transistors are powered from the M.V. power circuit.

The following is a description of the operation of this flip flop.

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In the absence of a pulse, transistor Q 3 is a conductor and Q 4 is blocked.

A positive pulse applied on the base of transistor Q 3 blocks it in its turn, producing on its collector a pulse of negative voltage which, being transmitted by one of condensers C 16 to C 19 on the base of Q 4, causes Q 4 to conduct.

The emitters of Q 3 and Q 4 are powered by the same resistance, so, when one conducts, the polarisation caused by the voltage drop in R 32 blocks the other.

This interdependence of the two transistors produces, at the moment of a change of state, a cumulative effect which reduces the rocking time.

The collector charge of Q 4 is constituted by one of potentiometers R 35 to R 39 in parallel on the Zener diode CR 11, the purpose of the latter being to clip the voltage pulse supplied by Q 4 at the moment when it becomes conducting, before sending it on to the integrating circuit.

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2-4 - INTEGRATION AND READING CIRCUIT :

In the clipped pulse supplied at the multivibrator outlet a certain energy is available of which a part is used to completely charge a capacity which is then discharged into a high-value integration condenser via a linearised integrating circuit containing diode CR 13 and transistor Q 5. The integration condenser thus receives a constant quantity of electricity at each hit counted.

The measuring device reads the charge voltage of the integration condenser.

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2-5 - POWER SUPPLY CIRCUITS :

Through the medium of transistor Q 2, transformer T 1 and the rectifying and stabilising circuits, the single set of batteries supplies the average stabilised at about 12 V for a current of 5 to 8 mA, as well as the high regulated voltages of - 500 v (50 μ A) and - 1 500 v (15 μ A) used by the photomultiplier.

The following summarises briefly the operation of the power supply :

The function of transistor Q 2 is to cut out, in the form of pulses, the current supplied by the batteries. These current pulses cross the primary winding of transformer T 1 and are kept at a regular rhythm by the reaction winding. The secondary windings feed pulses to rectifier circuits which, in the case of the H.V. secondary winding may also be voltage multipliers and which are followed by filtering and stabilisation circuits.

In the H.V. winding, a voltage of about 800 v is obtained by a first rectification and is then filtered before powering the 500 v corona-effect regulating tube VR 3. The high voltage of 1 500 v is obtained by a mono-alternation doubling of the voltage also collected on the H.V. winding of the transformer, which then supplies a voltage of 1 000 v connected in series with the 500 v. This 1 000 v voltage is also stabilised by a regulating tube VR 2.

The operating frequency of the converter is about 3 Kc/s

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2-6 - BATTERY-CHECKING CIRCUIT :

On the «Contrôle» position of the switch the micro-ampéremeter is connected as a voltage meter and enables the voltage at the terminals of the batteries to be checked. The measurement takes place under charge in order to avoid the influence of the internal resistance of the batteries. These are still usable when the needle goes inside the blue area.

2-7 - DIAL-LAMP CIRCUIT :

The dial lamp lights up in all switch positions except «Arrêt» (stop). This enables the operation of the equipment to be checked, the lamp being powered by the high voltage of the photomultiplier.

2-8 - SOUND ALARM CIRCUIT (buzzer) :

The sound alarm is controlled by a monostable flip flop producing signals that are rich in harmonics, whence the name «buzzer».

The recurrence frequency of the square pips emitted by the monostable flip flop is a direct function of the current passing into the measuring device of the SPP 2 NF and, consequently, of the counting rate.

This flip flop contains transistors Q 9 and Q 10. R 54 ensures a DC link between the collector of Q 9 and the base of Q 10. The time constant defining the duration of the pips is given by R 53 and C 36. Q 8 constitutes a gate controlling the release of the monostable. Q 11 is a common-emitted power stage necessary for correct operation of the sound unit. The C 34, R 50, C 35 assembly forms an integrating filter intended to weaken the fluctuations existing in the collector current of Q 5. Diode CR 18 polarised by the dividing bridge R 48-R 49 serves to clip the voltage appearing across R 50 on the Q 8 emitter in order to avoid saturation of the monostable when the control current reaches excessive values.

The following is a brief summary of the operation of the buzzer :

In the absence of signals from the detection probe, the collector current of Q 5 (situated on the «Integration» plate) is zero, and the four transistors Q 8, Q 9, Q 10 and Q 11 are non-conducting.

When counting signals operate the CR 13-Q 5 integrating circuit, transistor Q 5 delivers a current proportional to the rate of counting. Its collector current is divided between R 50, the base and collector circuits of Q 8 on the one hand and R 40, CR 20 on the other hand (the latter are on the battery container). With a current circulating in its base, Q 8 starts to conduct and its collector current charges C 36 across R 52 and R 53. When the voltage on the base of Q 9 reaches the V_{be} value, this causes Q 9 to conduct and the monostable Q 9-Q 10 to rock. While the pip in the monostable lasts, the three transistors Q 9, Q 10, Q 11 conduct and the membrane of the sound unit (LS 2) receives a pulse. At the same time, the collector current of Q 10, which sets up a voltage drop in R 56, leads a positive polarisation across R 51 on to the base of Q 8, which blocks this transistor.

At the end of the cycle, C 36 is charged with a positive on the collector side of Q 10. It rapidly discharges across R 56 and CR 19. When C 36 is discharged, Q 9, Q 10 and Q 11 are insulating and the original conditions are restored. Through its collector Q 8 again transmits the

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drive current across R 52 and R 53 on to C 36. The voltage appearing on the base of Q 9 will reach the V_{be} value more quickly according as the current supplied by Q 8 is high.

Consequently when the collector current of Q 5 increases which corresponds to arise in the counting rate, the pips in the monostable draw closer and the sound frequency given off by the unit increases.

It should be noted, moreover, that transistor Q 5 is inserted in a constant-current mounting, that is, for a given counting rate the collector current of Q 5 is stable whatever the collector charge may be, at least within the equipment's normal operating limits. The operational threshold of the monostable can therefore be moved by varying R 40, since in this way one changes the proportion of collector current flowing in R 40 + CR 20 and in R 50 + Q 8 + R 51 + R 56.

The threshold control button, that is, the button of potentiometer R 40 mounted in variable resistance, is situated on the battery container, and is accessible during use, as was seen above.

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1 - INITIAL STARTING :

At the time the equipment is first used it does not normally contain batteries. Proceed then as follows :

1-1 - FITTING OF BATTERIES :

Remove the lid of the battery container and place three 1,5 v torch batteries on the supports provided. Replace the lid.

1-2 - PRELIMINARY CHECK :

Place the button of the switch on the «contrôle» position ; the dial lamp of the measuring device should light up and the latter's needle should move into the blue area.

If the needle stays below this area, the batteries are defective or worn out and should be replaced.

2 - MEASUREMENTS OF RADIOACTIVITY :

Before any series of measurements is performed, the condition of the batteries must always be checked on the «contrôle» position.

Using the switch button, choose the 15 000 c/s position. If the intensity of the radioactive field in which the SPP 2 NF is placed is sufficient, the needle should move. After a few seconds it is stabilised, and then indicates the average number of real hits detected by the scintillator crystal in a second.

If the deviation of the needle is less than a third of the scale, go over to the 5 000 c/s position. The deviation will be greater and the reading accuracy increased. If the deviation is still less than one third of the scale, go over in succession to the positions 1 500, 500 and then 150 c/s until the needle is in the last two-third of the scale. It should be mentioned that the last sensitivity of 150 c/s is situated on the last two positions of the switch with two different integration time-constants permitting use of the device on the «Rapide» (fast) time constant during prospecting and on the «Lent» (slow) time constant during measurement.

3 - MEASUREMENT OF DOSAGE STRENGTHS IN mR/h :

As regards the direct gamma radiation of Cobalt 60, the correspondence between hits/second and dosage strength (mR/h) is directly indicated by the curve appearing on the reverse side of the equipment case.

When the average energy E of detected radiation is different from 1.25 MeV, the value indicated by the correspondence curve must be multiplied by the coefficient K read from the curve at the end of the present handbook Appendix 6.

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Determination of the coefficient K of course involves a knowledge of the average energy of the radiation in question.

4 - ADJUSTMENT FOR NORMAL SERVICE :

All settings of the electronic system of the SPP 2 NF scintillometer are carried out in the factory and should not be altered subsequently for a major repair. The calibration procedure is described later. Setting of the threshold of the sound alarm is performed by pressing the milled button of the potentiometer situated on the battery container beside the probe cable outlet. Identify, on the dial of the measuring device, the point where the sound alarm is activated, by bringing the control source close to the probe terminal and then moving it away.

5 - OPERATION

For use of the SPP 2 NF in the field, carrying equipment is provided composed of a set of cloth straps of adjustable length, fitter with snaps which hook on to the case, enabling it to be used at the waist.

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1 - MAINTENANCE FOR NORMAL SERVICE :

When the SPP 2 NF is used in regular service, the only maintenance necessary should be the replacement of the batteries when they become exhausted.

The operations for replacing the batteries are the same as those for the initial starting.

2 - STORAGE :

If an SPP 2 NF scintillometer is to remain unused for an extended period (a month or more) the power supply batteries should be removed from their housing and put in a dry cool place.

3 - REPAIR AND CALIBRATION :

Provided the equipment has not gone wrong the calibration need not be altered. However, when it has broken down and undergone a major repair, it may be necessary to perform re-calibration.

It should be noted that the calibration is more or less unaffected if a part is replaced in the power-supply and photomultiplier-crystal sub-assemblies.

The new calibration may then concern only the adjustment of the counting threshold.

If a part is changed on the electronic assembly (including the measuring device ; the ranges should be adjusted as follows :

Open the pistol probe by unscrewing the part covering the measuring device. Take out the electronic block and the photomultiplier scintillator block. Disconnect the terminals of the blue and white H.V wires on the electronic side. Leave in place the pulses from a generator into sleeve n° 7 of the socket formerly occupied by the photomultiplier (7th sleeve from the stop going counter clockwise, socket seen from the side opposite the welds).

These pulses must be positive with respect to earth, must have a width of 2 microseconds and an amplitude of about 2 volts.

With the switch on the 15 000 c/s position and the generator supplying 12 000 pulses per second, potentiometer R 35 is adjusted to obtain a reading of 12 000 c/s.

The adjustment potentiometers are situated on the «Integration» plate (see position on the wiring diagram)

The same procedure is carried out on all ranges, the measuring points being :

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When all the ranges have been adjusted, raise the threshold of the electronic assembly, which is about 50 mV on the base of transistor Q 3 of the flip flop . The counting threshold should then be adjusted with the use of a 137 Cs source. Beforehand, replace the photomultiplier scintillator assembly and the blue and white H.V wires. Replace the electronic assembly in the pistol body .

4 - ADJUSTMENT OF THE COUNTING THRESHOLD (see wiring diagram)

Bringing a 137 Cs close to the detector assembly, observe on the oscilloscope the pulses transmitted at the base of transistor Q 3 on the rocker plate. The pulses due to the photoelectric peak and shown as an excess are seen very clearly. These pulses correspond to the gamma rays of 660 KeV energy. Their level is adjusted with the oscilloscope by means of potentiometer R 30 (buzzer plate) in such a way that gamma rays of an energy less than 30 KeV are not counted. This level is given by the formula.

$$V_{mV} \text{ 137 Cs} = \frac{\text{Threshold of the electronics (about 50mV)} \times 660}{30}$$

It should be noted that potentiometer R 10 situated in the pistol (see diagram) allows coarse adjustment of this level by varying the gain of the photomultiplier. This adjustment is performed once and for all in the factory. It should only be altered if the photomultiplier is replaced. Once the level of the counting threshold is adjusted, the adjustments and calibrations are complete.

5 - QUICK CHECKING OF THE ADJUSTMENT :

A 137 Cs source (supplied with the equipment) gives a counting rate when it is applied, on the label side, against the front of the pistol.

The purpose is to permit a rapid check of the equipment. The method is the following :

With the source a long way from the equipment (at least 3 to 4 metres), the proper ambient movement is measured. The source is then applied against the front of the pistol. It is then checked that counting rate is approximately the sum of the proper movement and of the indication appearing on the source.

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










ACCESSORIES DIVERS
Casque
Echelle portative
Enregistreur

VARIOUS ACCESSORIES
Headset
Portable ladder
Recorder

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
CORRECTION FACTOR K
APPLICABLE TO CORRESPONDENCE
CURVE C/S ↔ mR/h
 (fixed on the instrument)

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